



Overview

Nylon planetary mill jar is a kind of high quality grinding jar which does not contain any metal impurity. It is the best option for milling all kind of metal and high level pigment. Nylon mill jar is often widely applied in many industries such as noble metal, electrical ceramics , magnetic materials, pigment coating, and so on.



Nylon planetary mill jar with agate balls



Nylon Planetary Mill Jar




Nylon planetary mill jar with zirconia balls

Feature

- 1.Pure white nylon material and elegant appearance
- 2.Non-metal material, no impurity in process of grinding
- 3.Professional design, good workmanship and beautiful outlook
- 4.High strength, good tenacity and tightness

Dimensions:



参数 (mm) 规格 (L)	D	D1	D2	D3	H	H1	H11	H2	b	δ	δ_1	R	O型圈规格
0.05	$\phi 63$	$\phi 43$	$\phi 48$		50	40	$+0.16$ 0	2.75	3.6	12	9	5	3.55 \times $\phi 47.5$
0.1	$\phi 68$	$\phi 48$	$\phi 53.5$		65	55	$+0.16$ 0	2.75	3.6	14	11	5	3.55 \times $\phi 53$
0.25	$\phi 88$	$\phi 68$	$\phi 73.5$		85	75	$+0.19$ 0	2.75	3.6	16	13	5	3.55 \times $\phi 73$
0.5	$\phi 108$	$\phi 88$	$\phi 93$		100	90	$+0.22$ 0	2.75	3.6	18	15	5	3.55 \times $\phi 92.5$
1	$\phi 133$	$\phi 109$	$\phi 115.5$		130	118	$+0.22$ 0	2.75	3.6	20	16	5	3.55 \times $\phi 115$
1.5	$\phi 133$	$\phi 109$	$\phi 115.5$		166	154	$+0.22$ 0	2.75	3.6	20	16	5	3.55 \times $\phi 115$
2	$\phi 158$	$\phi 134$	$\phi 140.5$	$\phi 145$	160	148	$+0.25$ 0	2.75	3.6	15	12	5	3.55 \times $\phi 140$
2.5	$\phi 158$	$\phi 134$	$\phi 140.5$	$\phi 145$	190	178	$+0.25$ 0	2.75	3.6	15	12	5	3.55 \times $\phi 140$

Parameter of Nylon Planetary Mill Jar

Model	Outside Diameter (mm)	Inner Diameter (mm)	Height (mm)
25ml	63	31	65
50ml	63	43	59
100ml	68	48	76
250ml	88	68	98
500ml	108	88	115
1L	133	109	142
1.5L	133	109	178
2L	158	134	172
2.5L	158	134	202
3L	158	138	210
4L	178	158	226
5L	218	190	217
10L	246	218	307
15L	273	243	372
20L	320	285	359
25L	320	285	424