





Overview

PU (Polyurethane) planetary mill jar has a superior performance in applications of electronics and battery materials due to its high wear resistance. It does not bring any harmful impurity into grinding materials and ensures high quality of electronic products.



PU Planetary Mill Jar



PU planetary mill jar with agate balls




PU planetary mill jar with zirconia balls

Feature

Tencan PU planetary mill jar has an integral structure with advantages of large circular arc, thick tank wall, easy opening, well-sealed, transparent tank wall, beautiful outlook, strong and durable. There are not any crack or split gap, no dead corner in process of grinding.

Dimensions



Parameters of PU Planetary Mill Jar

Model	Outside Diameter (mm)	Height (mm)
100ml	66	77
250ml	100	95
500ml	105	105
1000mL	135	140

