




Overview

Tempered nylon planetary mill jar is a kind of modified nylon jar by mold casting and integral forming, which has excellent performance in wear resistance, It is a kind of preferred mill jar for planetary ball mill duo to its durability.



Tempered Nylon Planetary Mill Jar



Tempered Nylon Planetary Mill Jar

Feature

1. Integral mold casting, good tightness, higher density than common nylon.
2. Good uniformity, no internal stress. press, more stable structure.
3. Tempered modification, high hardness, good stiffness, small friction
4. High hardness, good wear resistance.
5. Good stiffness, good impact resistance.
6. Lower friction coefficient, lower heat.



Tempered Nylon Planetary Mill Jar

Parameter of Tempered Nylon Planetary Mill Jar

Model	Outside Diameter (mm)	internal diameter (mm)	Height (mm)
250ml	98	73	90
500ml	109	84	120
1000ml	133	108	140
2000ml	159	145	175
3000ml	178	154	198
5000ml	219	188	218
10L	248	216	300

