



## Overview

Corundum, this name comes from India and it is a kind of mineral name. its hardness is only less than diamond and is mainly used for grinding all advanced materials, Corundum mill jar matched with zirconia balls is the best and classical solution for wear resistance. We'd like to remind that although corundum is very strong in hardness, it is easy to be broken. Therefore, you are requested to prevent from falling down to the ground or high strength impact to avoid break when using it. Any size of this kind of vacuum mill jar can be ordered and customized according to customer's requirements.



Corundum Planetary Mill Jar



Corundum Planetary Mill Jar




Corundum Planetary Mill Jar

## Feature

High temperature resistance, corrosion resistance, high strength, high hardness, wear resistance and etc.

## Dimensions



规格	D	D1	D2	D3	R	H	H1	H2	H3	O型圈规格
25ml	48	32	36.7	42	3	42	8	13	2	$\phi 36.5 \times 3.55$
50ml	55	40	44	49	3	50	8	13	2	$\phi 43.7 \times 3.55$
100ml	66	50	55	60	4	62	8	15	2	$\phi 54.5 \times 3.55$
250ml	86	70	75.5	80	4	80	8	16	2	$\phi 75 \times 3.55$
500ml	104	86	93	96	5	100	8	16	2	$\phi 92.5 \times 3.55$
1L	132	112	118.5	122	5	117	10	17	3	$\phi 118 \times 3.55$
1.5L	146	126	132.5	136	6	137	10	17	3	$\phi 132 \times 3.55$
2L	158	138	145.5	146	6	151	12	19	3	$\phi 145 \times 3.55$
2.5L	170	150	176	156	8	160	12	19	3	$\phi 175 \times 3.55$
3L	170	150	176	156	8	185	12	19	3	$\phi 175 \times 3.55$

## Parameter of Corundum Planetary Mill Jar

Model	Outside Diameter (mm)	Height (mm)
50ml	56	70
100ml	70	72
250ml	95	88
500ml	106	112
1000ml	133	131
1500ml	132	173
2000L	150	171