



## Overview

PU is abbreviation of polyurethane. PU is a kind of new high-tech and organic high polymer material. PU products are widely used in industries of electric, vehicle, construction, light industry, aerospace, textile, etc. [http://dict.youdao.com/javascript:void\(0\);](http://dict.youdao.com/javascript:void(0);) Pu roll mill jar has an excellent performance for grinding materials in these industries.



PU Roll Mill Jar



PU Roll Mill Jar



PU Roll Mill Jar

## Feature

1. Good wear-resistance, no scratch, no noise, long service life and low cost.
2. High mechanical strength, good elasticity, good recoverability.
3. Temperature resistance within  $-20^{\circ}\text{C}$  to  $120^{\circ}\text{C}$ .
4. Oil resistance and anti-aging, good weather resistance, long life long of 15 to 20 years.
5. No pollution, non-poison and no smell. [http://dict.youdao.com/javascript:void\(0\);](http://dict.youdao.com/javascript:void(0);)

## Parameter of PU Roll Mill Jar

Model	Outside Diameter (mm)	Internal Diameter (mm)	Height (mm)
100ml	68	53	69
250ml	98	74	90
500ml	106	87	104
1L	132	105	135
1.5L	132	116	161
3L	159	139	223
5L	192	168	251
10L	243	219	286

